PATHFAST™ TB LAM Ag Control

[English]

Issue Date: 01/06/2024

Revision: ver.2

Intended use

PATHFAST TB LAM Ag Control is used for quality control of PATHFAST TB LAM Ag chemiluminescent immunoassay with the PATHFAST automated analyser. PATHFAST TB LAM Ag Control is to be used by laboratory technicians, nurses or physicians, in a hospital or clinical laboratory setting.

Package composition of provided materials

Control 1 (CTL-1)	2.0 mL x 3 bottles (liquid, CMIT/MIT* 0.0015 - < 0.01%)
Control 2 (CTL-2)	2.0 mL x 3 bottles (liquid, CMIT/MIT* 0.0015 - < 0.01%)
CONTROL DATA SHEET	1 sheet
Instructions for use	1 sheet

*CMIT/MIT: Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (3:1)

Classification according to Regulation (EC) 1272/2008: Skin sensitisation (Category 1A). For hazard and precautionary statement, see "Chemical Hazard & Precautionary statements".

Precautions and warnings

- Dispose of all measured reagents and materials according to the standard disposal method. For example, autoclave at 121 °C for 20 minutes. Follow general precautions, and handle all components as if capable of transmitting infectious agents.
- This product should not be used for any other purpose.
- When any serious incident occurs in relation to the product, report it to the manufacturer and the competent authority in the area which the user and/or the patient are located.

Preparation before use

Ready to use.

Storage and expiration

- Store at 2 - 8 °C.
- The expiration date is listed on each bottle and package.
- Avoid contamination and exposure to direct sunlight.
- Avoid water damage during storage.
- After opening, controls can be used until the expiration date.
- Do not use controls beyond the indicated expiration date.

Quality Control assay (QC assay)

Please refer to "Preparation and procedure" in the instructions for use of PATHFAST TB LAM Ag.

- Use CTL-1 or CTL-2 as a sample.
- Place the reagent cartridge in the cartridge rack, and then dispense approximately 100 µL of the control into a sample well and load onto PATHFAST.
- Carry out sample assay in accordance with PATHFAST operator's manual.

Target values and ranges

The target values are traceable to the in-house standard. The target values and ranges were determined and evaluated by PHC Corporation. They were obtained using PATHFAST TB LAM Ag assay reagents and an analyser available at the time. The target values and ranges are listed on the "CONTROL DATA SHEET" packaged with PATHFAST TB LAM Ag Control.

* PATHFAST: JP Registered Trademark No. 5982733

* PATHFAST: US Registered Trademark No. 3074207

Chemical Hazard & Precautionary statements

Warning!

May cause an allergic skin reaction.

Avoid breathing mist. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of soap and water.

Contains Mixture of 5-chloro-2-methyl-4-isothiazolin-3-one and 2-methyl-4-isothiazolin-3-one (3:1) [EC No. 911-418-6, CAS No. 55965-84-9].

Supplier:

PHC Corporation

1460-6 Aza-Mitodai, Mito, Tako-machi,

Katori-gun, Chiba, 289-2247 Japan

+31 76 543 3833 (c/o PHC Europe B.V.)

Summary of safety and performance is available from:

European Database on Medical Devices (EUDAMED).

Contact for technical assistance

www.pathfast.eu/contact

Symbols

PHC Corporation uses the symbols and signs on the last page in addition to those listed in EN ISO 15223-1:2021 (Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements).

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