PHC Corporation	DEVICE MASTER FILE PATHFAST Presepsin Control Product No. PF0201C	Document No.	YM-R02-7020CD02
		Revision No.	4

Chapter 2 INFORMATION TO BE SUPPLIED BY THE MANUFACTURER 2.2 Instructions for use
--

1)	English	p. 1
2)	German	p. 1
3)	Italian	p. 2
4)	Spanish	p. 2
5)	Romanian	p. 3
6)	Czech	p. 4
7)	Dutch	p. 4
8)	French	p. 5
9)	Greek	p. 5
10)	Hungarian	p. 6
11)	Lithuanian	p. 7
12)	Polish	p. 7
13)	Slovenian	p. 8
14)	Japanese	p. 8

PATHFAST™ Presepsin Control

[English]

Issue Date: 01/06/2024

Revision: ver.2

* PATHFAST: JP Registered Trademark No. 5982733

* PATHFAST: US Registered Trademark No. 3074207

Intended use

PATHFAST Presepsin Control is used for quality control of PATHFAST Presepsin chemiluminescent immunoassay on PATHFAST automated analyser. PATHFAST Presepsin Control is used by a laboratory technician, nurse or physician, in hospital including the emergency room, doctor's office and clinical laboratory. PATHFAST Presepsin Control is a device for near patient testing (NPT).

Package composition of provided materials

Control 1 (CTL-1)	For 1.0 mL x 2 vials (lyophilized)
Control 2 (CTL-2)	For 1.0 mL x 2 vials (lyophilized)
Control diluent	1.0 mL x 4 bottles (liquid, Na azide < 0.1%)
Instructions for use	1 sheet

Precautions and warnings

- Azide can react with the copper and lead used in some plumbing systems to form explosive salts. When disposing of azide-containing materials, they should be flushed away with large volumes of water.
- Dispose of all measured reagents and materials according to the standard disposal method; for example, autoclave at 121 °C for 20 minutes. Follow general precautions, and handle all components as if they are capable of transmitting infectious agents.
- This product should not be used for any other purpose.
- If any serious incident occurs in relation to the product, report it to the manufacturer and to competent authorities in the area in which the user and/or the patient is located.

Preparation before use

Transfer whole volume of one bottle Control diluent totally into one vial control. Stand closed for 5 minutes to reconstitute. Mix carefully. Avoid foam.

Transfer aliquots of the freshly reconstituted controls into appropriate tubes for storage. Store the aliquots immediately at -20 °C or lower. Perform only one control procedure per aliquot.

Note:

Use the same lot of CTL-1 or CTL-2 and Control diluent when reconstituting a control. Never mix different lots of CTL-1 or CTL-2 and Control diluent.

Storage and expiration

- Store at 2 - 8 °C.
- The expiration date is listed on each bottle and package.
- Avoid contamination and exposure to direct sunlight.
- Avoid water damage during storage.
- Do not use reagents and controls beyond the indicated expiration date.

Stability of the reconstituted control

At 2 - 8 °C	28 days
At -20 °C or lower	6 months
After thawing	Use only once

Quality Control assay (QC assay)

Please refer to the "Preparation and procedure" in instructions for use of PATHFAST Presepsin.

- Use a reconstituted CTL-1 or CTL-2 as sample.
- Place the reagent cartridge in the cartridge rack, and then dispense approximately 100 µL of the control into a sample well and load onto PATHFAST.
- Carry out sample assay in accordance with PATHFAST operator's manual.

Target values and ranges

The target values are traceable to in-house standard material assigned by amino acid analysis. The target values and ranges were determined and evaluated by PHC Corporation. They were obtained using PATHFAST Presepsin assay reagents and analyser available at the time. The target values and ranges are listed on "CONTROL DATA SHEET" packaged in PATHFAST Presepsin.

Summary of safety and performance is available from:
European Database on Medical Devices (EUDAMED).

Contact for technical assistance

www.pathfast.eu/contact

Symbols

PHC Corporation uses the following symbols and signs in addition to those listed in EN ISO 15223-1:2021 (Medical devices - Symbols to be used with information to be supplied by the manufacturer - Part 1: General requirements).

This symbol means "Device for near patient testing".

(Symbols for self-testing and near-patient testing under the IVD Regulation 2017/746/EU. MedTech Europe. Dec. 13, 2018)

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PHC Europe B.V.
Axon Lab AG