


# 3D Cell Culture in Gravity TRAP™

**KEYWORD** 1) Spheroid 2) 3D culture 3) Hanging Drop 4) Label free assay

**SUMMARY** Spheroids formed by the method of hanging drop were moved to Gravity TRAP™ ULA plate, and quantified by Cell<sup>3</sup>iMager neo.


## Materials and Methods

Cell Line: HepG2 cells (RIKEN BRC)  
 NIH3T3 cells (RIKEN BRC)  
 Medium: DMEM (Nacalai tesque)  
 Plate: GravityTRAP™ ULA plate (InSphero)  
 Seeding cell density: 500, 1000cells / well  
 Culture days: 5 days after spheroid formation by hanging drop  
 Imaging methods: Bright-field, 4800dpi  
 Bracket Focus (stacked)




## Results and Conclusions

- It was possible to quantify the spheroids formed by hanging drop for growth rate estimation.



**Pseudo Volume Sum (a.u.)** Growth rates of HepG2 spheroids were higher than NIH3T3 spheroids.



## SCREEN Holdings Co., Ltd.

KYOTO(Head office) / Tenjinkita 1-1, Teranouchi-agaru 4-chome, Horikawa-dori,  
 Kamigyo-ku, Kyoto 602-8585, Japan

**Life Science Business Development and Sales Division**

KYOTO(Rakusai)

Furukawa-cho 322, Hazukashi, Fushimiku, Kyoto 612-8486, Japan

Phone : +81-75-931-7824 / Fax : +81-75-931-7826

TOKYO

7th Floor, Yamatane Bldg., 2-21 Etchujima 1-chome, Koto-ku, Tokyo 135-0044, Japan

Phone : +81-3-4334-7977 / Fax : +81-3-4334-7978

<http://www.screen.co.jp/eng>

International Callers;

Call: +1 847 870 7400 X 2423 (or) +1 847 910 3374